



Farbe	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Liner	1	2	2	2	4	1	1	2	2	1	3	3	4	1	5	5
Titan Opaquer	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Core Dentin	CD2	CD6	CD6	CD3	CD5	CD6	CD6	CD3 CD5 ¹⁾	CD5	CD2	CD4	CD5	CD5	CD4	CD5	CD5
Opaque Dentin	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Dentin	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Enamel	LIGHT YELLOW	YELLOW	YELLOW	DARK ORANGE	DARK ORANGE	LIGHT YELLOW	YELLOW	ORANGE	DARK ORANGE	YELLOW	YELLOW	ORANGE	DARK ORANGE	LIGHT YELLOW	YELLOW	ORANGE
Enamel²⁾	57	57	58	59	59	57	58	58	59	58	59	59	59	59	59	59



1) CD3 | CD5 im Verhältnis 50:50

2) Bei der Verwendung der Core Dentine kommen die Enamelmassen 57-59 zum Einsatz.